

# MCS<sup>®</sup>+ 8150

The gold standard in  
red cell collection

INTEGRATED DEVICES



# The gold standard in red cell collection



With the proven reliability and mobility of the MCS<sup>®</sup>+ device, collecting multiple blood components from a single donor to meet collection goals has never been easier. With the ability to run multiple protocols on the same device, the MCS+ system gives operators the flexibility to collect the components in the greatest demand.

# Features and benefits

## Provides proven, reliable technology

- An industry leader with more than two million red cells collected
- Simple to operate
- Increases donor management capability – uses fewer donors to meet your collection needs

## Versatility: Multiple protocols to meet collection needs

- 2 units of red blood cells
- 2 units of red blood cells (leukocytes reduced)
- 1 unit red blood cells and a single 400–550 mL unit of plasma
- 1 unit red blood cells and a jumbo plasma that can be split into 2 smaller units

## Mobility: Reach more donors

- Compact, lightweight, mobile weighing only 56 lbs
- Stackable transport case for maximum protection in the mobile environment

## Revision L.1 Software

- New female donation requirement:
  - 5'3" and > 150 lbs
  - Opportunity to convert additional type O and B female whole blood and RBCP donors to 2RBC for incremental red cell unit
- Additional operational and software enhancements addressed with new software

## Implementation support to ensure program success!

Expert support for program implementation and extensive post-installation monitoring delivers enhanced staff productivity to meet collection and operational goals



## To contact Customer Service...

Phone: (800) 537-2802 ■ Fax: (800) 860-1512

Email: [CustomerServiceNA@haemonetics.com](mailto:CustomerServiceNA@haemonetics.com)

| Description  | List No.    | Unit / Case |
|--|-------------|-------------|
| MCS®+ 8150 device  | 08150-110-E | 1           |
| Two unit red cells, leukocytes reduced, disposable                   | 832F-00     | 8           |
| Two unit red cells disposable  | 832-00      | 20          |
| Red cell and plasma disposable                                       | 822-00      | 20          |
| Red cell and plasma disposable (plasma can be split post collection) | 822-2P      | 8           |
| Transport case   | 09012-00    | 1           |

### Device Specifications

| Approx. dimensions | in operation | closed      |
|--------------------|--------------|-------------|
| Height             | 26.5 inches  | 17 inches   |
| Width              | 21.5 inches  | 21.5 inches |
| Depth              | 21.5 inches  | 12 inches   |



**RxOnly**

© 2006, 2008-2009, 2012 Haemonetics Corporation.  
Haemonetics and MCS are trademarks or registered trademarks  
of Haemonetics Corporation in the USA, other countries, or both.  
01.2012, Germany. COL-FS-000059(AC)